



**Steve Mosey**, Co-owner Mobility Technologies, LLC 262 Fort Wayne Avenue Richmond, IN 47374

lifetrike.com 765.983.8711 info@lifetrike.com

[2] Gill Sans font



Helvetica font



Steve Mosey, Co-owner Mobility Technologies, LLC 262 Fort Wayne Avenue Richmond, IN 47374

lifetrike.com 765.983.8711 info@lifetrike.com

Helvetica font



Steve Mosey, Co-owner Mobility Technologies, LLC 262 Fort Wayne Avenue Richmond, IN 47374

lifetrike.com 765.983.8711 info@lifetrike.com



Mobility Technologies, LLC 262 Fort Wayne Avenue Richmond, Indiana 47374







LifeTrikeEnvelope.indd 1





- Life Trike Sample Documents

Mobility Technologies, LLC 262 Fort Wayne Avenue • Richmond, Indiana 47374-2328 www.lifetrike.com • 1-765-983-8785







Mobility Technologies, LLC 262 Fort Wayne Avenue • Richmond, Indiana 47374-2328 www.lifetrike.com • 1-765-983-8785

# EXPERIENCE THE ULTIMATE FREEDOM IN MOBILITY WHILE DRIVING THE ...





This inno ative vehicle changes lives!





Mobility Technologies, LLC

262 Fort Wayne Avenue • Richmond, IN 47374

info@lifetrike.com

lifetrike.com

Call: 1.765.983.8711





#### STANDARD FEATURES

TWIN ENGINE, Low center of gravity, Low carbon footprint, Emergency ramp/gate release, Twin turn signal bulbs for redundancy on left and right, Twin fuel tanks, Twin batteries, Will run on one engine (reduces horsepower by half), Adjustable handlebars, Tight turning radius (6' 6"), Twin caliper front brakes, Reverse gear, 12 volt power socket

### SPEC SHEETS

- Turning radius: 6 feet 6 inches
- · Length: 72 Inches
- · Width: 67 Inches
- Height: 43 Inches (47 inches to top of mirrors)
- GVW: 544 lbs without occupant. Lifetrike.com
- · Occupant weight limit: 250 lbs.
- · Wheel chair weight limit 256 lbs.
- Total Weight limit 506 lbs
- · Maximum wheel chair size 18 inch. (26 inch across wheels)
- Ground clearance: 3-5 Inches
- Tire pressure front and rear: 30-36 psi
- Engine size 50cc x 2
- Engine style: 4 Stroke
- Open road MPG: 52 City 45
- · Required fuel: 92 octane gasoline

Get Out and Ride







Experience the Ultimate Freedom in Mobility while driving the most Innovative Vehicle available today. Introducing the



STABLE and SECURE ride.

EASY ON EASY OFF ELECTRIC RAMP. CRUISE up to 30 mph.

NO RIDER transfer Necessary.

ACCOMODATES Powered and WHEELCHAIR TRANSPORTS.

ADJUSTABLE STEERING controls fit all ages, heights, and ABILITIES.

DRIVER comfort ASSURED.

MAXIMIZE LOCAL TRAVEL WITH EXTENDED RANGE.

**CLICK ON OUR VIDEOS AND PUT YOURSELF** BEHIND THE WHEEL AND FEEL THE WIND IN YOUR FACE

LifeTrike.com

### EXPERIENCE A NEW LEVEL OF MOBILITY COMFORT WITH LifeTrike.

DESIGNED AND ENGINEERED BY MOBILITY TECHNOLOGIES, LLC.

View our videos and watch how the **LifeTrike** offers you increased mobility and range of travel. Visit

## LifeTrike.com



**OUR STORY** 







### How can I get financial assistance?

call us direct at (765) 983-8711.

o I buy the LifeTrike?

Check with your Dealer for financing options and promotional discounts.

our website for the Dealer in your area, or

**WLifeTrike** 

TLY ASKED QUESTIONS

### How do I get my *LifeTrike* fixed?

All local moped and scooter Dealers can repair your vehicle. Our technicians are available to offer Factory Support when necessary. Custom engineered repair parts are stocked and shipped from our Richmond, IN manufacturing facility.

### Is my *LifeTrike* street legal?

Yes. All *LifeTrikes* come with factory-equipped lights, turn signals, horn, and brake lights to comply with all municipal enforcement requirements and regulations. Most states require some form of registration or licensing.



### SERVICE CENTERS

- Repairs performed by Local Moped and Scooter Dealers.
- Repair parts are available through Dealer network.
- Repair Manuals & Parts Lists available upon request.
- Customized Adaptations and Engineering Designs provided by Mobility Technologies, LLC.









sidewalks toward our destination we noticed a gentleman about a half a block ahead of us who was in a power wheelchair trying to negotiate his way around some trash cans that had been placed too far out on the sidewalk. Well the next thing to happen was lot good thing. He dropped one of his wheels of of the curb and toppled over into the treet. A horrifying sight to say the least. We rulhed up to help him up off of the sirc and and get him on his way. Fortunate who was uninjured

but I realized at that ment hat this incident

could have turner out much worse.

In March of 2008 my wife and I were vaca-

As many may know who have been to Key

West you will generally find that walking

almost anywhere in the town is possible.

to take a walk from our hotel to a place on

one of the main streets for some breakfast.

tioning with some friends in Key West Florida.

On one sunny day during our stay we decided

While walking along one of the narrow

As we returned home on the plane
I couldn't get the image out of my hear of
how vulnerable the person was or the toe
Clearly he needed a safer my to trave

I began to think altout the schoters that we rented to sightsee the day after the accident and it occurred to me that if somehow you could harry the liberating sense of freed me schoter offers with the limits of a wheelchair then that person could be off of the sicewalk and on the street where he cants to be in the first place.

In act, after I returned home I began to notice things about power chairs that mad previously not noticed; the dangers and frustrations of trying to navigate busy city streets. This was worse I thought than what I saw in Key West. Something had to be done about this. I shared my story with my brother, George Mosey, who was equally disturbed. Through a collaborative effort, the LifeTrike was born!

~ Steve Mosey, President

Life it was strikes the right balance between comfort and safety to provide riders a thrilling and enjoyable experience. Extend the capability and range of your power chair or wheelchair to enhance your lifestyle by choosing when and where you'll travel.

The LifeTrike offers an increased level of personal freedom, and is the only affordable, limited mobility product on the market today. \*\*

Made in the USA.

This innovative vehicle changes lives!





